

Product name: **Anti- Bovine Albumin**

Product code: **AB8031**

Immunogen: Bovine serum albumin

Host species: Chicken

Antibody type: Polyclonal

Serum fraction: IgY fraction

Conjugate: _____

Formulation: Phosphate buffered saline (PBS) containing sodium azide* (10 mM) as the preservative

Storage: **At 4-8 °C.** If keep in freezer, avoid repeated thawing and freezing

**Concentration/Titer: Protein: 15 mg/ml

**Sodium azide is very toxic if swallowed. Contact with skin, wash immediately with plenty of water.*

** Working dilution in each experiment should be determined by titration assay.

This product is for laboratory use only

www.razibitech.com

info@razibitech.com

