

Product name: **Anti-Goat IgG FITC conjugate**  
Product code: **AF8038**

Immunogen: Goat IgG whole molecule

Host species: Rabbit

Antibody type: Polyclonal

Serum fraction: IgG fraction

Conjugate: Fluorescein isothiocyanate (FITC)

Formulation: Phosphate buffered saline (PBS) containing carrier protein and sodium azide\* (10 mM) as the preservative

Storage: At 4-8 °C, **Do not freeze**

\*\*Titer: Indirect fluorescent assay: 1:50 – 1:120

*\*Sodium azide is very toxic if swallowed. Contact with skin, wash immediately with plenty of water.*

*\*\* Working dilution in each experiment should be determined by titration assay. Due to differences in assay systems, the titer that we have defined may not reflect the user's actual working dilution.*

**This product is for laboratory use only**

www.razibiotech.com



info@razibitech.com