

Product name: **Anti-Human Immunoglobulins**

Product code: **AH7078**

Immunogen: Human Immunoglobulins ( $\gamma$ ,  $\alpha$ ,  $\mu$ )

Host species: Goat

Antibody type: Polyclonal

Serum fraction: IgG fraction

Conjugate: \_\_\_\_\_

Formulation: Phosphate buffered saline (PBS) containing sodium azide\* (10 mM) as the preservative

Storage: **At 4-8 °C.** If keep in freezer, avoid repeated thawing and freezing

\*\*Concentration/Titer: Protein: 15 mg/ml

*\*Sodium azide is very toxic if swallowed. Contact with skin, wash immediately with plenty of water.*

\*\* Working dilution in each experiment should be determined by titration assay.

**This product is for laboratory use only**

www.razibitech.com

info@razibitech.com

