

Product name: **Anti-Human Immunoglobulins**

Product code: **AH7181**

Immunogen: Human Immunoglobulins (H+L)

Host species: Goat

Antibody type: Polyclonal

Serum fraction: IgG fraction

Conjugate: _____

Formulation: Phosphate buffered saline (PBS) containing sodium azide* (10 mM) as the preservative

Storage: **At 4-8 °C.** If keep in freezer, avoid repeated thawing and freezing

**Concentration/Titer: Protein: 15 mg/ml

**Sodium azide is very toxic if swallowed. Contact with skin, wash immediately with plenty of water.*

** Working dilution in each experiment should be determined by titration assay.

This product is for laboratory use only

www.razibitech.com

info@razibitech.com

