

Product name: **Anti-Human Immunoglobulins**  
**( $\gamma$ ,  $\alpha$ ,  $\mu$ ) HRP conjugate**

Product code: **AP7077**

Immunogen: Human Immunoglobulins  
( $\gamma$ ,  $\alpha$ ,  $\mu$ )

Host species: Goat

Antibody type: Polyclonal

Serum fraction: IgG fraction

Conjugate: Horseradish peroxidase (HRP)

Formulation: Phosphate buffered saline (PBS) containing carrier protein and preservative (ProClin 300)

Storage: At 4-8 °C, **Do not freeze**

\*Titer: ELISA: 1:10000 – 1:40000  
Western blot: 1:4000 – 1:15000  
Dot blot: ND

\* Working dilution in each experiment should be determined by titration assay. *Due to differences in assay systems, the titer that we have defined may not reflect the user's actual working dilution.*

**This product is for laboratory use only**

www.razibiotech.com

info@razibitech.com

