

Product name: **Anti-Human Immunoglobulins (H+L) HRP conjugate**

Product code: **AP7088**

Immunogen: Human Immunoglobulins (H+L)

Host species: Goat

Antibody type: Polyclonal

Serum fraction: IgG fraction

Conjugate: Horseradish peroxidase (HRP)

Formulation: Phosphate buffered saline (PBS) containing carrier protein and preservative (ProClin 300)

Storage: At 4-8 °C, **Do not freeze**

*Titer: ELISA: 1:15000 – 1:40000
Western blot: 1:4000 – 1:20000
Dot blot: ND

** Working dilution in each experiment should be determined by titration assay. Due to differences in assay systems, the titer that we have defined may not reflect the user's actual working dilution.*

This product is for laboratory use only

www.razibiotech.com

info@razibitech.com

