

Product name: **Anti-Human IgA HRP conjugate**
Product code: **AP7180**

Immunogen: Human IgA₁ and IgA₂ whole molecules

Host species: Goat

Antibody type: Polyclonal

Serum fraction: IgG fraction

Conjugate: Horseradish peroxidase (HRP)

Formulation: Phosphate buffered saline (PBS) containing carrier protein and preservative (ProClin 300)

Storage: At 4-8 °C, **Do not freeze**

*Titer: ELISA: 1:10000 – 1:40000
Western blot: 1:4000 – 1:20000
Dot blot: ND

* Working dilution in each experiment should be determined by titration assay. *Due to differences in assay systems, the titer that we have defined may not reflect the user's actual working dilution.*

This product is for laboratory use only

www.razibitech.com

info@razibitech.com

