

Product name: **Anti-Fish IgM HRP conjugate**  
Product code: **AP8005**

Immunogen: Fish IgM whole molecule

Host species: Rabbit

Antibody type: Polyclonal

Serum fraction: IgG fraction

Conjugate: Horseradish peroxidase (HRP)

Formulation: Phosphate buffered saline (PBS) containing carrier protein and preservative (ProClin 300)

Storage: At 4-8 °C, **Do not freeze**

\*Titer: ELISA: 1:15000 – 1:50000  
Western blot: 1:4000 – 1:15000  
Dot blot: ND

*\* Working dilution in each experiment should be determined by titration assay. Due to differences in assay systems, the titer that we have defined may not reflect the user's actual working dilution.*

**This product is for laboratory use only**

www.razibitech.com

info@razibitech.com

