

Product name: **Anti- Bovine IgG**
Product code: **AB7083**

Immunogen: Bovine IgG whole molecule

Host species: Rabbit

Antibody type: Polyclonal

Serum fraction: IgG fraction

Conjugate: _____

Formulation: Phosphate buffered saline (PBS) containing sodium azide* (10 mM) as the preservative

Storage: **At 4-8 °C**. If keep in freezer, avoid repeated thawing and freezing

**Concentration/Titer: Protein: 10 mg/ml

**Sodium azide is very toxic if swallowed. Contact with skin, wash immediately with plenty of water.*

**** Working dilution in each experiment should be determined by titration assay.**

This product is for laboratory use only

www.razibiotech.com

info@razibiotech.com

