

Product name: **Anti-Mouse IgG**
Product code: **AM7088**

Immunogen: Mouse IgG whole molecule

Host species: Goat

Antibody type: Polyclonal

Serum fraction: IgG fraction

Conjugate: _____

Formulation: Phosphate buffered saline (PBS) containing sodium azide* (10 mM) as the preservative

Storage: **At 4-8 °C.** If keep in freezer, avoid repeated thawing and freezing

**Concentration/Titer: Protein: 15 mg/ml

**Sodium azide is very toxic if swallowed. Contact with skin, wash immediately with plenty of water.*

**** Working dilution in each experiment should be determined by titration assay.**

This product is for laboratory use only

www.razibitech.com

info@razibitech.com

